

Integrated
Innovative

**SMART
CALIBRATION
FACILITY**

For
HbA1c, CRP, RF, ASO
in cost effective
way

h of Class in Biochemistry...

Features

- Direct access to more than 200 tests through test keys
- Maintenance Free peristaltic pump
- Unique circuitry for long lamp life
- Triple cuvette system
- Real time graph with online reading
- Multi standard calibration & memory
- Levey jennings and standard deviation graphs
- Displays actual temperature on screen
- Automatic calibration against standard
- Sample carry over prevention facility
- Robust system with built in stabilizer
- Patient report with PID & name
- User Friendly software
- USB interface

Technical Specifications

OPTICAL SYSTEM

Linear measurement range : 0.000 to 3.000 absorbance units (A)

Photometric accuracy : $\pm 2\%$ or 0.007 whichever is higher, From 0 to 1.5 A, $\pm 3\%$ From 1.5 A to 3.0 A

Optical measurement : Photodiode

Light Source : Tungsten halogen

ANALYSIS MODES

Biochemistry & Turbidimetry : Absorbance, End Point,

Fixed Time, Kinetic, Differential & Ratio

Coagulation : PT, INR & APTT

FILTERS

Type of filter : Hard coated narrow band interference

Wave length : 340, 405, 510, 546, 578, 630 nm & two optional

FLOW CELL

Sipping volume : upto 1000 μ l

Triple cuvette system : 18 μ l flow cell

10 mm square cuvette

Round tube for coagulation

POWER

Wattage : 100 VA

Voltage : 115 - 230 Volts $\pm 10\%$, 50 - 60 Hz

INCUBATOR

15 Positions at 37°C

MEMORY

Open tests : 250 tests

Results : 2500 patient results

QC : 30 QC results/test

DISPLAY

Human machine interface : TOUCH PANEL/KEYPAD
5" Color Graphics LCD

PRINTER

On board thermal 28 columns, maintenance free

STORAGE CONDITIONS


Temperature : From - 10°C to + 60°C

Relative humidity : Up to 85%


Enclosure : ABS Fire retardant

Dimensions in mm (lxbxh) : 380 x 305 x 175 (Apr.)


Weight : 5 kgs. (Apr.)




Online Graph



LJ/SD QC Chart



Incubator



Printer

Due to continuous development, specifications are subject to change without prior notice.

invitronixTM
quality precision service

Gala No. 61/62, 3rd Floor, Building B, Gami Industrial Park, TTC MIDC Industrial Area, Pawne,
Navi Mumbai - 400705. Email : info@invitronix.com, Web : www.invitronix.com